

# Chemistry 2.1

# Quantitative Analysis

Titration Tasks

# Titration experiment #1- Standard solution

**Aim:** Prepare a standard solution of  $\sim 0.1 \text{ molL}^{-1}$  of  $\text{Na}_2\text{CO}_3$

**Type:** Procedure

**Method:**

- 1) Calculate the mass of  $\text{Na}_2\text{CO}_3(\text{s})$  require to produce 250 mL of  $0.1 \text{ molL}^{-1}$   $\text{Na}_2\text{CO}_3$  solution.
- 2) Weigh the calculated amount of  $\text{Na}_2\text{CO}_3(\text{s})$  into the weighing cup using electronic balance.
- 3) Record the exact mass
- 4) Half fill your volumetric flask with distilled water
- 5) Transfer the solid into the flask using a filter funnel
- 6) Rinse the cup as well as the funnel into the flask at least twice to ensure all substance is transferred
- 7) Fill the volumetric flask to the mark

## Titration experiment #2- Standardising HCl

**Aim:** To find the concentration of unknown HCl

**Type:** Quantitative

### Method:

- 1) Fill the burette with standard  $\text{Na}_2\text{CO}_3$  solution
- 2) Pipette 20 mL of unknown HCl to a flask.
- 3) Add two to three drops of indicator to the flask.
- 4) Add  $\text{Na}_2\text{CO}_3$  into the flask until end point reached
- 5) Record the titre, repeat process until 3 titres are within 0.2mL range.

Trial	#1	#2	#3	#4
<i>Titre (mL)</i>				

- 6) Using the results, calculate the concentration of unknown HCl

# Titration experiment #3- Finding unknown NaOH

**Aim:** To find the concentration of unknown NaOH

**Type:** Quantitative

**Method:**

- 1) Fill the burette with standard HCl solution
- 2) Pipette 20 mL of unknown NaOH to a flask.
- 3) Add two to three drops of indicator to the flask.
- 4) Add HCl into the flask until end point reached
- 5) Record the titre, repeat process until 3 titres are within 0.2mL range.

Trial	#1	#2	#3	#4
<i>Titre (mL)</i>				

- 6) Using the results, calculate the concentration of unknown NaOH

# Titration experiment #4- Finding unknown $\text{H}_2\text{SO}_4$

**Aim:** To find the concentration of unknown  $\text{H}_2\text{SO}_4$

**Type:** Quantitative

## Method:

- 1) Fill the burette with standard NaOH solution
- 2) Pipette 20 mL of unknown  $\text{H}_2\text{SO}_4$  to a flask.
- 3) Add two to three drops of indicator to the flask.
- 4) Add NaOH into the flask until end point reached
- 5) Record the titre, repeat process until 3 titres are within 0.2mL range.

Trial	#1	#2	#3	#4
<i>Titre (mL)</i>				

- 6) Using the results, calculate the concentration of unknown  $\text{H}_2\text{SO}_4$